

## FACS Staining Protocol

### Materials

Reagents of known titration ( $\lambda$  per 45 $\lambda$  rxn volume)

Cells, washed, adjusted to 2.5-5 x 10<sup>7</sup>/ml or 5-10 x 10<sup>5</sup>/well

3911 Falcon flexible plastic plate

Staining media (SM): PBS/3% FCS/0.04% NaN<sub>3</sub>. keep on ice.

Centrifuge that can spin plates

Ice

Falcon 2052 tubes

### Procedure:

1. Make up staining mixes, and second step e.g.

Reagent	$\lambda$ /rxn	# rxns	total vol
RS3.1-bi	0.3	10	3.3
a-mu-fl 1	10	10	
SM	25 total	10	250-10-3.3=237

Should make up one more rxn than needed. Note, remember to include single color and unstained controls. The single color controls should not get PI. Make one unstained tube with PI and one without.

2. Spin at 4° C about 13K xg for 20-30 min to deaggreagate.

3. Wash cells into SM. Adjust cell count to 2.5-5 \* 10<sup>7</sup>/ml. Can also aliquot approx. number of cells directly to plate, then wash with 2 x 100 $\lambda$  of SM, reconstituting in about 20 $\lambda$

4. If preparing cells separately, add 25 $\lambda$  of rxn mix to plate first. Then add 20 $\lambda$  of cells to each appropriate well, mixing and changing tip as you go. By using a crossing pattern, can add same stains to various cell preps or various cell preps to same stains.

5. If preparing cells in plate, add stains directly after cells are washed and reconstituted.

6. On ice for 20', covered with foil.

7. Add 100 $\lambda$ /well of SM, using multichannel pipettor. Spin 1500 rpm x 3' at 4° C.

8. Suck off sup using 9" pasteur hooked to gentle suction. Tilt plate at 45° angle to visualize pellet. Enter the tip of pipette into well at a 45° angle along the edge. Do not try to suck dry.

9. Tap plate gently to dislodge pellets. Add 100 $\lambda$  SM. Spin. Suck. Wash 1x more again.

10. If using streptavidin or second step, add 50 $\lambda$ /well, mixing and changing tips as you go. This should have been prepared in #1 above. Otherwise, skip to step 13.

11. On ice for 15', covered.

12. Wash as in steps 7-9.

13. Final resuspension in 100 $\lambda$  of SM. Transfer to labelled (I use the 96wp grid A-H, 1-12 for labels) 2052 tubes which contain 200 $\lambda$  of SM. If using PI, then add at appropriate concentration (e.g. 2ug/ml final- dilute 1:500 of our stock) to the SM before dispensing to

Mark Shlomchik  
September 14, 2004

the 2052 tubes. Transfer cells to the 2052 tubes. Keep on ice and in the dark. Analyze on  
FACS